

Bed Stick



- Provides a support pole for the side of the bed
- Ideal for patients requiring assistance while turning, repositioning during transfers
- Slides between mattress and bed base
- Can be used on either or both sides of the bed
- Easily installed and removed
- Powder coated steel

Specifications:

- Height: 510mm
- Reach: 350mm

Lithgow Mobility
25 Main Street
Lithgow NSW 2790
02 6352 2752 / 0428 090 993

Bathurst Mobility
112 Keppel Street
Bathurst NSW 2795
02 6332 2370 / 0484 113 223