

# Bed Wedge



Large folding bed wedge

An Icare Bed Wedge can be used to elevate the upper body to alleviate acid reflux, respiratory problems, and back or neck pain.

Icare Bed Wedge can also be used to increase circulation and reduce conditions such as varicose veins in the legs.

## **Ventilated Base Fabric**

Vant fabric is a three-dimensional polypropylene material that has high air permeability and is anti-bacterial. Using this in the base greatly increases air movement.

### **Specifications:**

- Zip off washable cover
- Three in one design
- Width: 78cm
- Length: 68cm
- Depth: 21cm (thick end)

The Icare bed wedge can also be placed under one end of the mattress for alleviating head or legs

Lithgow Mobility  
25 Main Street  
Lithgow NSW 2790  
02 6352 2752 / 0428 090 993

Bathurst Mobility  
112 Keppel Street  
Bathurst NSW 2795  
02 6332 2370 / 0484 113 223