

# Removable Bed Rail



- Provides a support rail on the side of the bed
- Ideal for patients requiring assistance while turning, repositioning or during transfers
- Easily installed and removed
- Slides between mattress and bed base
- Powder coated steel

## Specifications:

- Width: 600mm
- Height: 500mm
- Reach: 590m

Lithgow Mobility  
25 Main Street  
Lithgow NSW 2790  
02 6352 2752 / 0428 090 993

Bathurst Mobility  
112 Keppel Street  
Bathurst NSW 2795  
02 6332 2370 / 0484 113 223