

## Swift Self Propelled



The Days Swift Self-Propelled Wheelchair with Handbrakes has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair. Self-Propelled in three different sizes.\*

### Features:

- Light frame 7.8kg
- Folding backrest (for storage or transport)
- Quick release wheels (for transport)
- Fixed arms
- Swing away and removable leg rests
- Solid tyres (no maintenance required)
- Handbrakes for added safety\*
- Extended brake handle (ease of brake application)
- Maximum User Weight: 115kg

Lithgow Mobility  
25 Main Street  
Lithgow NSW 2790  
02 6352 2752 / 0428 090 993

Bathurst Mobility  
112 Keppel Street  
Bathurst NSW 2795  
02 6332 2370 / 0484 113 223

<b>Specifications</b>	<b>18 x 16"</b>
Overall Height:	930mm
Overall Width:	660mm
Overall Depth:	1090mm
Seat Width:	450mm
Seat Depth:	425mm
Seat Height at Front:	500mm
Seat Height at Rear:	475mm
Seat Angle Adjustment:	N/A
Seat Belt:	No
Upholstery:	Black, nylon, padded, fire retardant
Push Handles:	Pair of wheelchair grips
Push Handle Height:	1080mm
Headrest:	N/A
Backrest:	Folding
Backrest Height from Seat:	400mm
Backrest Angle Adjustment:	Fixed 90°
Armrests:	Polyurethane foam, fixed, plastic side/skirt guards
Armrest Height from Seat:	220mm
Frame:	Aluminium
Frame Colour:	Silver
Frame Weight:	7.8kg
Type:	Self-Propelled
Legrest Length:	420 – 540mm
Knee Angle Adjustment:	N/A
Leg/Footrest:	Swing away, 120mm adjustable length (tool required), removable
Leg/Footrest Angle Adjustment:	Fixed 70°
Footplate:	215 x 150mm, lift up
Front Fork:	Plastic
Front Wheels:	200mm (8"), solid polyurethane
Rear Wheels:	610mm (24"), solid polyurethane, quick release
Brakes:	Hand brake / Attendant brake, direct to rear wheel/parking brake
Net Weight:	12.5kg
Maximum User Weight:	115kg
Warranty:	12 months
ARTG ID:	174573
Tested to Australian Standards:	AS/NZS 3695.1:2011 Wheelchair standard

Lithgow Mobility  
25 Main Street  
Lithgow NSW 2790  
02 6352 2752 / 0428 090 993

Bathurst Mobility  
112 Keppel Street  
Bathurst NSW 2795  
02 6332 2370 / 0484 113 223