

## Leg Lifter



- This webbing strap with stiffened stem enables the user to move a stiff or immobile leg, when using the bed, wheelchair, footstool, etc.
- The upper loop fits around the hand and wrist, and the lower stiffened loop slides over the foot to give control when lifting the leg.
- Length 864mm.
- Weight 78g.

Lithgow Mobility  
25 Main Street  
Lithgow NSW 2790  
02 6352 2752 / 0428 090 993

Bathurst Mobility  
112 Keppel Street  
Bathurst NSW 2795  
02 6332 2370 / 0484 113 223